KS KEIR
SURGICAL

Instructions For Use
Instructions d'utilisation

Instructions for Use

Initial Use

- All instruments which are being used for the first time must be thoroughly cleaned, rinsed and then sterilized before using.

Cleaning

- Soiled instruments should be soaked in a detergent or moistened using a pre-clean foam/gel. Instruments should not soak for long periods in tap water nor should blood be allowed to dry on the instruments.
- Whenever possible, instruments that can be dismantled should be taken apart for cleaning
- An ultrasonic bath is appropriate for hard to clean instruments.
- Machine or manual cleaning is acceptable. Care should be taken to assure that detergents are compatible with instrument materials.
- Following chemical disinfection, instruments should be rinsed using distilled water.
- Hinged instruments should be opened before cleaning.
- Considerations for the weight and bulk of instruments in a basket should be made during the cleaning process to ensure that cleaning fluids can circulate freely.
- Deposits remaining after the cleaning process should be completely removed before sterilization.
- Instruments with many moving parts or specific articulation should be lubricated with a steam soluble lubricant (milk).
- Instruments must be carefully dried after cleaning and rinsing is complete.

Testing for Correct Function

- Instruments should be carefully and thoroughly checked before each use.
- After each cleaning, instruments should be visibly free from protein-based residue and other soil.
- Care should be taken so that any cracks, breakage or corrosion is noticed.
- Instruments with flecks should be separated and taken out for special treatment.
- All moving parts as well as working tips and blades should be checked with special care.
- Instruments which are in any way defective or damaged should be immediately sorted out from the others and sent to Keir Surgical for repair or replacement.

Sterilization

- All stainless steel instrumentation from Keir Surgical may be steam sterilized based upon the following criteria:
 - 270°F Pre-Vacuum Steam Cycle (132°C)
 - Exposure Temperature: 270° – 275°F (132° – 135°C)
 - Exposure Time: 4 minutes
- Instruments with ratchets or self-retaining mechanisms should be sterilized in the open position.
- Please follow the autoclave manufacturer's operating instructions. All machines used must comply with all current quality standards.
- Sterilization packing material or containers must comply with current standards.

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