**KATENA  
SURGICAL INSTRUMENTS  
HANDLING INSTRUCTIONS  
Please read before use!**

**English**

**Caution**

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

All persons using this device should be knowledgeable in the use and handling of surgical instruments, accessories and related equipment.

**Initial Use of New Instruments**

Every instrument must be cleaned and sterilized before it is used for the first time.

**Inspection and Functional Check**

It is very important to carefully examine each surgical instrument for breaks, cracks or malfunctions before use. It is especially essential to check areas such as blades, points, ends, stops and snaps as well as all moveable parts.

Do not use damaged instruments. Never attempt to make repairs yourself. Service and repairs should be referred to trained qualified persons only. Refer questions about repair to the manufacturer.

**Cleaning and Maintenance**

Every surgical instrument should be disinfected and thoroughly cleaned after each use. Inspect and test each instrument carefully. Sterilize all instruments before surgery. A good cleaning and maintenance procedure will ensure correct function and extend the useful life of the instrument.

Special attention should be paid to slots, stops, ends, hollow tubes and other highly inaccessible areas. Check insulation, cables and connectors for cuts, voids, cracks, tears, abrasions, etc.

**Do not use damaged instruments.**

If an instrument can be disassembled, it should be cleaned in the disassembled state.

(Note: Keep all parts together and protected from being misplaced, do not intermix with other or similar devices of other manufacturers.)

Cleaning and rinsing must take place immediately after each use for best effect. Failure to clean promptly may result in adherent particles or dried secretions that may resist cleaning and complicate or resist future sterilization. Instruments must be completely cleaned and rinsed of all foreign matter.

Use warm water and a commercially available instrument pre-soak or cleaning agent. Enzymatic cleaners should be used to remove protein deposits. If used, follow the enzymatic cleaner's instructions and rinse thoroughly.

- Do not use corrosive cleaning agents (i.e. bleach). Cleaning solutions and rinses at or near a neutral pH (7.0) are best.
- Do not use abrasive cleaners.
- Only a soft bristle brush should be used.
- Flush all inside channels and brush with a small diameter bottle or tube brush.
- Ultrasonic cleaners can be used. Check and retighten any fittings that may have vibrated loose.
- Rinse thoroughly with distilled water.
- Prepare for storage and/or sterilization.

After cleaning and rinsing, dry instruments completely and carefully with compressed air (including inside channels and highly inaccessible areas).

Note: After cleaning and before sterilization, it is strongly recommended that all moving parts, lock boxes, joints and catches be lubricated with a physiologically safe lubricant.

**Storage and Sterilization**

Instruments should be stored in a clean, dry, moisture free area. Micro surgical instruments should be stored individually in their shipping carton or in a protective tray with partitions. Always keep a protective cover on delicate tips when the instrument is not in use.

Katena instruments are reusable and meet ANSI/AAMI standards for sterilization. Steam autoclave sterilization is recommended. Thoroughly clean instruments of all debris, tissue and foreign matter prior to sterilization.

Follow the sterilizer manufacturer's instructions for operation and loading of steam autoclaves. There must be direct steam exposure to all surfaces of the instruments being sterilized including the internal surface and tubes channels.

Allow instrument to air cool to room temperature before use.

**Standard Sterilization Methods**

Steam autoclave sterilization is recommended. Standard autoclave cycle:

Steam sterilize at 270°F (132°C) for 5 minutes. Other time and steam temperature cycles may also be used. However, user must validate any deviation from the recommended time and temperature.

(Note: Contact the manufacturer of your steam autoclave to confirm appropriate temperatures and sterilization times.)

Autoclave temperatures should not exceed 280°F (137°C); handles, insulation or other non-metallic parts may be damaged.

**Handling**

All surgical instruments should be handled with the greatest care when being transported, cleaned, treated, sterilized and stored. This is especially true for blades, fine points and other sensitive areas.

Surgical instruments corrode and their functions are impaired if they come into contact with aggressive materials. The instruments should not be exposed to acids or other aggressive cleaning agents.

**Warranty**

All Katena surgical instruments are unconditionally guaranteed against manufacturing defects. Any instrument found to have a manufacturing defect will be repaired or replaced free of charge.

Katena's surgical instruments are designed and manufactured to meet the highest quality standards. We cannot accept any liability for failure of products which have been modified in any way from their originals or for any failure due to misuse or application which is not in accordance with the designer's intentions.

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