

# Finger Splints

## Finger Splints Stack Type



**Indications:** Mallet finger, fingertip and nail bed injuries

**Skintone and Clear Thermoplastic:** Used to fix the distal finger joint in the extension position. The hood over the phalanx media should preferably be fixed using the Velcro® strap. The proximal interphalangeal joint remains unobstructed. The side walls of the splint are perforated for ventilation.

Stack Type		-skintone or clear thermoplastic-		
Skintone Item No.	Clear Item No.	Size	Length (cm)	Finger width (cm)
10-2225/10	10-2228/10	00	4.5	3.5
10-2225/11	10-2228/11	0	4.5	4.0
10-2220	10-2228/01	1	5.5	4.5
10-2221	10-2228/02	2	5.3	5.0
10-2222	10-2228/03	3	6.0	5.5
10-2223	10-2228/04	4	6.2	6.0
10-2224	10-2228/05	5	6.0	6.5
10-2224/55	10-2228/55	5 <sup>1</sup> / <sub>2</sub>	6.0	7.0
10-2225	10-2228/06	6	6.5	7.5
10-2225/07	10-2228/07	7	6.2	8.0

## Finger Splints Split Splints



**Indications:** Mallet finger, finger tip and nail bed injuries

**Adjustable:** The adjustable, modified finger splints are split right along the dorsal proximal section. Consequently, when the fixing strap is put into position, the splint fits this area of the finger exactly. As the swelling of the finger subsides, the splint may be adjusted easily.

Finger Splints		-skintone-	
Item No.	Size	Length (cm)	Finger width (cm)
10-2232/10	00	4.5	4.5
10-2232/11	0	4.5	4.5
10-2232/01	1	5.5	5.5
10-2232/02	2	5.3	5.3
10-2232/03	3	6.0	6.0
10-2232/04	4	6.2	6.2
10-2232/05	5	6.0	6.0
10-2232/55	5 <sup>1</sup> / <sub>2</sub>	6.0	6.0
10-2232/06	6	6.5	6.5
10-2232/07	7	6.2	6.2

## "STACTILE" Finger Splints Stack Type



**Indications:** Mallet finger, finger tip and nail bed injuries

**Stactile, Skintone and Clear Thermoplastic:** In view of the anatomy, the STACTIP® Finger Splints have a distal, volar window over the tactile elevation. Thanks to this window, it is possible for the patient to touch and to carry out certain sensitive tasks with "finger tip sensation".

Finger Splints		-skintone or clear thermoplastic-		
Skin Tone Item No.	Clear Item No.	Size	Length (cm)	Finger width (cm)
10-2208/01	10-2209/01	1	5.4	4.5
10-2208/02	10-2209/02	2	5.3	5.0
10-2208/03	10-2209/03	3	5.5	5.5
10-2208/04	10-2209/04	4	6.0	6.0
10-2208/05	10-2209/05	5	5.8	6.5
10-2208/55	10-2209/55	5 <sup>1</sup> / <sub>2</sub>	6.0	7.0
10-2208/06	10-2209/06	6	6.4	7.5
10-2208/07	10-2209/07	7	6.4	8.0

## Boutonnière Finger Splints



**Indications:** Boutonnière Deformity

**Boutonnière, Skintone and Clear Thermoplastic:** The Stack type Finger Splints have been specially designed for fingers with "buttonhole deformity" (Boutonnière Deformity).

Finger Splints		-skintone or clear thermoplastic-		
Skintone Item No.	Clear Item No.	Size	Length (cm)	Finger width (cm)
10-2218/01	10-2219/01	1	5.0	4.5
10-2218/02	10-2219/02	2	5.0	5.0
10-2218/03	10-2219/03	3	5.0	5.5
10-2218/04	10-2219/04	4	5.0	6.5
10-2218/05	10-2219/05	5	5.4	7.0
10-2218/06	10-2219/06	6	5.4	7.5