

Safely Retrieve Needles and Blades

MAGNATREVE® *Magnetic Sweep*

***Reduces The Risk of Deadly Needle
Sticks to Healthcare Professionals***

Minimize Employee Exposure to
Bloodborne Pathogens

Safely Retrieve Contaminated Needles
and Blades

Use as Part of Your Exposure
Control Plan / Engineering Controls

No moving parts to break, no wheels to fall off.
Twelve powerful magnetic heads (floor sweep)
are in direct contact with the floor allowing
hands-free collection of even the smallest
needles and ferrous metal objects.



Hands-Free Sharps Collection

MAGNATREVE®

*Use as Part of Your Exposure
Control Plan / Engineering Controls*

Magnetic Sweeps are Available in Three Sizes



100 TEENSY SWEEP

8"length x 2 1/2"head
(20.3 cm x 6.3 cm)



103 HAND HELD

12"length x 7 1/4"head
(30.5 cm x 18.4 cm)



105 FLOOR SWEEP

36"length x 12 1/2"head
(91.4 cm x 31.7 cm)

**Safely Retrieve
Needles and Blades**



Free with every order.

MAGNATREVE® forceps
safely remove sharps from
magnetic heads.

100% money back guarantee
Ask about volume purchase discounts

ADEPT+MED®
INTERNATIONAL

Working to Make Surgery Safer

Adept-Med International, Inc.
665 Pleasant Valley Road
Diamond Springs, CA 95619

*Prevent room-to-room contamination,
keep a MAGNATREVE® in every suite.*

**KS KEIR
SURGICAL**

604.261.9596 | 1.800.663.4525
www.keirsurgical.com

© Registered Trademark of Adept-Med International, Inc.

The MAGNATREVE® is Produced Exclusively in the USA

2/09 24M