



Chroma**PRO**TM

LED Light Source

POWERFUL, SOLID-STATE LED

- 50-Watt LED Illumination
- 50,000-Hour Lamp Life – Eliminates Bulb Changes
- 6,000°K Daylight-White Color Temperature
- 4 Port Rotating Turret – BFW/ ACMI, WOLF, OLYMPUS, STORZ
- Half the Energy to Produce Twice the Output of Halogen



Intelligently bright. Very bright.

TM

Chroma**PRO**[™] LED Light Source



- 50-Watt LED Illumination
 - Patent Pending Collection Optics and Cutting Edge LED Technology
 - 50,000 Hour Solid-State LED Limits Output Decline to 15% over Lifetime
 - Analog Dimmer Controls Output from 10% to 100%
 - Universal Power Supply
 - Thermal Protection Via On-Board Diagnostics
 - Small Footprint
 - Optimizes Light Transmission Through Cable
 - CE, IEC 60601-1, UL 60601-1, CSA C22.2 No 601.1, IEC 60601-1-2
- Chroma**PRO**[™] SYSTEM (Shown)**
INCLUDES XtreamBeam**[™] Fiber-Optic Headlight**
- A Fiber Optic Headlight Ergonomically Designed for Comfort, Low Profile and Durable Construction
 - Clean, Even, Adjustable Spot .5" to 4.25" @ 16" Working Distance (1.2cm – 10.8cm @ 40 cm Working Distance)
 - 4.75mm Diameter Cable with 9,000 Strands Enables Use with 30% Damage

PRODUCT CODES

BFW 5500 Chroma**PRO**[™] LED Light Source

BFW 9850 Xtream**Beam**[™] Fiberoptic Classic Headlight w/ 7' Cable

BFW 9851 Xtream**Beam**[™] Fiberoptic Sport Headlight w/ 7' Cable

BFW 9855 Xtream**Beam**[™] Fiberoptic Classic Headlight w/ 10' Cable

BFW 9856 Xtream**Beam**[™] Fiberoptic Sport Headlight w/ 10' Cable

BFW 9886 Mobile Stand for Chroma**PRO**[™]

KS KEIR
SURGICAL

604.261.9596 | 1.800.663.4525
www.keirsurgical.com



Intelligently bright. Very bright.