



CoAXTM High-Definition Headlight Imaging System


- True Coaxial Viewing
- 1080i or 720p Output
- HDMI or S-Video, Composite, Component, SDI Output
- CoAX^{TM-r} Model Offers Optional USB Port for Recording SD Still Shots and Video



Intelligently bright. Very bright.

CoAX™

High-Definition
Headlight Imaging System



Advanced System Features Include:

- **1080i or 720p Output** – Use with any HDMI monitor or recording device. Also compatible with S-Video, Composite or Component SDI Output
- **Incremental Digital Zoom** - 4x to 8x magnification of the surgical field
- **Incremental Digital Sharpness Settings** – Enhance actual viewing work window to the most diminutive variance in available light while producing high definition color and sharpness
- **One-Touch White Balance** – Easy, simple and accurate in seconds
- **Auto Exposure Setting** – User can isolate and enhance a portion of the viewing area with increased definition and depth of field when needed or for use in cases with minimal incision
- **True Coaxial Viewing** – What the surgeon sees is always on camera – Simply point and shoot without adjustment during a case
- **CoAX™-r Model Offers Optional USB Port for Recording SD Still Shots and Video** – Use Thumb Drive for immediate “Take-Away” imaging of case
- **Innovative Solid-State Plasma Lighting** – System illuminated by our ChromaLUME™ solid-state plasma light source – providing in excess of 10,000 hours of sunlight equivalent (5,800° K) illumination without lamp changes

Product Codes

BFW 7930HD – (NTSC) CoAX™ High-Definition Headlight Imaging System includes: CoAX™ Headlight/Camera Unit, 10' fiber optic cable, CoAX™ Control Unit with HDMI, S-Video, BNC and Composite output, ChromaLUME™ Solid State Plasma Light Source, Mobile Stand, Video Cables

BFW 7935HD – Same as above in PAL format

BFW 7930HDR – (NTSC) CoAX™ High-Definition Headlight Imaging System as listed above with built in USB port capable of capturing video and stills (as pictured to the left)

BFW 7935HDR – Same as above in PAL format

KS KEIR
SURGICAL

604.261.9596 | 1.800.663.4525
www.keirsurgical.com



Intelligently bright. Very bright.