



TOR-MARK™

by kmt

All-in-one Toric Marking Device

Product Code: TOR-MARK-10
Tor-Mark by KMT, Box of 5

Invented by Canadian ophthalmologists, the patent pending Tor-Mark™ is a single-use device for creating precise angular markings on the cornea for use in toric IOL implantation and L.R.I. procedures.

- ✓ Bubble level for accuracy
- ✓ Easily set marking angles for one-step marking technique
- ✓ Marking blades create clear, accurate marks
- ✓ Sterile / single use

Distributed by:

KS KEIR
SURGICAL

T: 1.800.663.4525 | 1.604.261.9596
E: info@keirsurgical.com | W: keirsurgical.com

Manufactured by:

CONKIN
SURGICAL INSTRUMENTS LTD.

T: 604.638.2757 | W: conkinsurgical.com