ENDO
KEIR SURGICAL



Featuring:
PMC™ Brushes

Flexible Endoscopy

2023 Edition

PRODUCTS FOR FLEXIBLE ENDOSCOPY

At Keir Surgical, we continually strive to meet the needs and challenges of Canadian hospitals and clinics by providing high quality products at the best possible price. Our goal is Customer Satisfaction.

FLEXIBLE ENDOSCOPY

A complete product line for endoscope care:

- Enzymatic bedside pre-clean kits, cleaning brushes, tip protectors, transport trays, drying and cleaning verification equipment.

Devices and accessories for GI and Bronchoscopy:

- Bite Blocks, Biopsy Forceps, Cytology Brushes, Needles, Snares, Retrievers & Polyp Traps.

All our Endoscopy products – Guaranteed to be LATEX-FREE.

GENERAL INFORMATION

Updated: February 2023

Ordering

Order by fax 604.261.9549, email orders@keirsurgical.com or EDI via Global Healthcare Exchange (GHX). For questions, please call 1.800.663.4525.

Shipment & Terms

Our products are shipped FOB for orders above $1000, unless otherwise specified. FOB does not apply to remote shipping locations. Payment terms are net 30 days. Due to fluctuating costs, Keir Surgical reserves the right to change prices without notice. Prices in effect at time of order prevail.

Returned Goods

To return merchandise, a Returned Goods Authorization number must first be obtained from our customer service department at 1.800.663.4525. All items returned must be either unopened or sterilized. Keir Surgical will not accept used, unsterilized items for return. All returned items are subject to inspection by Keir Surgical and/or the product manufacturer prior to issue of a credit or replacement. Return of unopened, unmarked packages will only be accepted within 30 days of invoice date. All returns are subject to a 20% restocking fee unless Keir Surgical and/or the Manufacturer deem product as defective in workmanship or a distribution error.

Product Information & Warranty

All our products are latex-free. The manufacturer warrants these products are free from defects and are manufactured with extremely high-quality workmanship. The storage conditions and the use of the products are factors beyond the control of Keir Surgical which may directly affect it and the results obtained from it. The user assumes responsibility for any consequences, which may arise from the final use of the products.

All reusable devices are warranted to the original purchaser by the manufacturer to be free from defects in materials and workmanship. All disposable devices are warranted to the original purchaser by the manufacturer to be free from defects in material and workmanship for one use only. The liability of products distributed by Keir Surgical does not extend to any misuse of these products.

We welcome your ideas and comments.

If you have a suggestion for a new device or an improvement of an existing device, please feel free to communicate your ideas to us.

Disclaimer

All brush recommendations made on the part of Keir Surgical are suggestions only. Due to the variety of instruments available on the market, end-users are in the best position to select the most suitable brush for their needs and shall assume responsibility for their choices.

*Please note that any images or pictures contained in our catalogue are for product demonstration purposes only and are in no way intended to be instructional or provide authority on acceptable use of products.

TABLE OF CONTENTS

ENDOSCOPE CLEANING & CARE

Inspection, Testing & Drying 6
Flexible Endoscope Brushes 9
Bedside Pre-Clean Supplies 14
Scope Cleaning & Tip Protection 15
Wipes & Cleaning Pads 16
Scope Transport Containers 17
Countdown™ Scope Tags 21

GI DEVICES & ACCESSORIES

Polyp Traps 24
Bite Blocks 25
Capsule Endoscopy 26
Retrieval Devices 27
Retrieval Forceps - Specialty 28
Polypectomy Snares 29
Injection Needles 30
Endoscopic Hemoclip 31
Cytology & Microbiology Brushes 33
Balloon Dilation System - GI 34
Balloon Dilators 35
Achalasia Dilators 37
Biliary Dilators 38
Biliary Stents 39
Pancreatic Flexi-Stents 40
Catheters 42

PRODUCT INDEX 43

Featuring:

PMC™ Brushes

KS KEIR SURGICAL

ENDOSCOPE CLEANING & CARE

INSPECTION, TESTING & DRYING 6

Inspection Scopes 6

Instrument Channel Dryer 7

FLEXIBLE ENDOSCOPE BRUSHES 9

PMC Brushes: Valve Brushes & Kit Combos 10

Flexible Endoscope Brushes 11

Valve Management 13

BEDSIDE PRE-CLEAN SUPPLIES 14

SCOPE CLEANING & TIP PROTECTION 15

WIPES & CLEANING PADS 16

SCOPE TRANSPORT CONTAINERS 17

Transport Carts 18

Transport Trays 19

Transport Bags 20

COUNTDOWN™ SCOPE TAGS 21

Featuring:

PMC™ Brushes

Keir Surgical Endoscope Care

KS COMPLETE CARE

Do complicated Endoscope Repairs cause you headaches?
Let us simplify that process.
We are pleased to offer flexible endoscope repairs, all in one convenient program.

OEM CORRECT™



WHAT WE OFFER:

- **Fast turnaround:** Repairs completed at our state of the art facility in Vaughan, ON
- **High Quality:** ISO 13485 certified & our team of technicians have over 250 years of combined OEM Correct™ endoscope repair experience!
- **Cost savings:** Significant savings compared to OEM repairs.

Comprehensive **Save Time & Money** **Trusted**

For more information, contact us at:
705.201.1188 | repairs@keirsurgical.com


INSPECTION SCOPES

Medical grade flexible borescopes


Video Scope

What's hiding inside your instruments?


Damage



Debris



Discolouration



Moisture



Visualization - Illumination - Magnification - Documentation







CODE	DESCRIPTION	COMPATIBLE CHANNEL
CD-USB	Dual Inspection Control Module, USB output	
CD-HDMI	Dual Inspection Control Module, HDMI output	
CB2-60	Inspection Borescope, OD 1.83 MM, 60 CM working length	> 2.0 mm
CB2-110	Inspection Borescope, OD 1.83 MM, 110 CM working length	> 2.0 mm
CB2-200	Inspection Borescope, OD 1.83 MM, 200 CM working length	> 2.0 mm
CB1-110	Inspection Borescope, OD 1.06 MM, 110 CM working length	> 1.1 mm

Note - software included for picture and video capture.

INSTRUMENT CHANNEL DRYER

Features

HEPA Filtered Air .01 micron with filter replacement alert.

Dry two scopes simultaneously with two independent pumps & timers.

Line pressure channel sensors immediately shut off Airtime pump and alert to potentially clogged channels.

Variety of mounting options - Wall, Pole Mount & Counter Stand.

Adjustable Dry times -10-minute default.

SPECS:
- Control Unit Dimensions: 10" H x 6" W x 4" D
- Control Unit Weight: 4.4 lbs
- Electrical Information: 100 - 240 VAC, 50 – 60 Hz
- UL & CE Marked
- DoE Energy Efficiency Level VI
- EC 60601-1 Medical Grade
- Other Features: Adjustable Timer Settings

airtime

CYGNUS MEDICAL



VFR100



Adapters & Tubing

Maximum flow rate reduces drying times.

Reduce risk of cross-contamination - tubing changed monthly and adapter sets after each use.

Compatible with Olympus, Pentax, Fuji, Storz & ACMI.

CODE	DESCRIPTION	PKG QTY
VFR100	Airtime™ Flexible Scope Dryer	1
DF100X	Tubing -Universal Directional Flow	20
SDT01	Tubing -Small Diameter Scope – Single line	5
SDT02	Tubing - Small Diameter Scope – Dual line	5
ASO3	Adapter Set– OLYMPUS- Single use	120
ASF2	Adapter Set – w/aux. H20 jet - FUJI - Single use	120
ASP2	Adapter set – w/aux. H20 jet - PENTAX - Single use	120
LL01-01	Luer-Lok™ Adapter & Cap (individual bags)	100
LL01-25	Luer-Lok™ Adapter & Cap (4 bags of 25)	100
FLT-Y30	HEPA Filter Assembly set - 0.010 mcg filtration	2
FLT-Y40ST	HEPA Filter Straight Assembly Set, 0.010 mcg	2


Non-sterile

INSTRUMENT CHANNEL DRYER

Mounting Kits

airtime

CYGNUS MEDICAL



IV Pole Mounting Kit

Cradle bar holds two scopes
Height reduction bars available
for non-adjustable poles



Wall Rail System

Dryer/Two scope set up
Height reduction system
3’ 6’’ wall space needed



Counter Stand & Bracket



Wall Mount Bracket

Mounts one Drying unit

CODE	DESCRIPTION	PKG QTY
IV500	Airtime™ IV Pole Mounting Kit - 1 Scope Cradle Bar 2 Height reduction bars, 1 Mounting Bracket & IV Pole	1
WRC100	Airtime™ Wall Mount Bracket - 12” rail, 2 wall mount anchors 1 universal wall rail mount bracket, 1 rail clamp	1
WRC500	Airtime™ Wall Rail System - Four 3-ft rails (white end caps) 2 scope cradles, 4 height reduction bars, 12 wall mount anchors 5 universal wall rail mount brackets & rail clamps	1
CS9179	Counter Stand with Universal Bracket	1

Non-sterile

SELECTING APPROPRIATE FLEXIBLE SCOPE BRUSHES

with: PMC™ Brushes

TAKE THE GUESSWORK OUT OF CHOOSING BRUSHES

Ask about our Complimentary Brush Consultation

We can help you choose the most effective brushes for your Flexible Endoscope models from the variety of sizes and styles we offer.

TIPS FOR CHOOSING A BRUSH:

1

CHANNELS
Identify scope channels/ports requiring brushing for each model.

2

TECHNIQUE
Consider brushing method required by MIFU - number of swipes; push or pull through.

3

APPROPRIATE SIZE
CSA standards focus on brush size for proper channel cleaning and damage prevention.



Confirm channel specs of scope make/model are compatible with brush min/max channel diameter & length. Refer to the corresponding column in PMC product charts.

See following example...

SCOPE TYPE	WORKING CHANNEL
Colonoscopes Standard/Therapeutic	Diam. 2.8 - 3.2 mm Length 135 - 165 cm
Gastroscopes Slim/Standard/therapeutic	Diam. 2.0 - 4.2 mm Length 110 - 140 cm
Bronchoscopes Intubation/Standard/Therapeutic	Daim. 1.2 - 2.8 mm Length 60cm
Cystoscopes Standard	Diam. 2.0 mm Length 65 cm
Ureteroscopes	Diam. 1.2 - 1.5 mm Length 67 cm


CODE	SHEATH COLOR	BRUSH DIA.	COMPATIBLE CHANNEL	LENGTH	SCOPE TYPE
MICROCHANNEL SCOPES					
KIT-0202-100	White	2.0 mm	1.0-1.7 mm	100 cm	Microchannel
SMALL DIAMETER SCOPES					
KIT-0303-120	Purple	3.0 mm	2.0-2.2 mm	120 cm	Bronch/Cystos
KIT-0303-230	Purple	3.0 mm	2.0-2.2 mm	230 cm	Slim/Ped. Gastros
GI SCOPES					
KIT-0505-230	Blue	5.0 mm	2.0-4.5 mm	230 cm	GI Scopes
KIT-0505-250	Blue	5.0 mm	> 2.8 mm	250 cm	Enteroscopes

Note - all about kits also include a 5/11mm valve & port cleaning brush.

INFO@KEIRSURGICAL.COM | 1.800.663.4525

PMC BRUSHES: VALVE BRUSHES & KIT COMBOS

Flexible Endoscope, Single use



Single-Ended and Double-Ended brushes for efficient cleaning of residues.
Channel friendly plastic & coated brush tip to avoid channel scratches.
Colored catheter, quick identification of different brush sizes.
Easy to use, individually packed kits.
Includes channel brush & combination port/valve brush.

VALVE CLEANING BRUSHES

CODE	DESCRIPTION	BRUSH DIAMETER	LENGTH	PKG QTY
VCB-10	Single-Ended	10 mm	7 cm	100
VCB-0511-15	Double-Ended	5 mm & 11 mm	15 cm	100

SINGLE ENDED BRUSH KITS

CODE	SHEATH COLOR	BRUSH DIAMETER	COMPATIBLE CHANNEL	LENGTH	SCOPE TYPE	PKG QTY
MICROCHANNEL SCOPES						
KIT-0002-100	White	2.0 mm	1.0-1.7 mm	100 cm	Microchannel	50
SMALL DIAMETER SCOPES						
KIT-0003-120	Purple	3.0 mm	2.0-2.2 mm	120 cm	Bronch/Cystos	50
GI SCOPES						
KIT-0005-230	Blue	5.0 mm	2.0-4.5 mm	230 cm	GI Scopes	50

DOUBLE ENDED BRUSH KITS

CODE	SHEATH COLOR	BRUSH DIAMETER	COMPATIBLE CHANNEL	LENGTH	SCOPE TYPE	PKG QTY
MICROCHANNEL SCOPES						
KIT-0202-100	White	2.0 mm	1.0-1.7 mm	100 cm	Microchannel	50
SMALL DIAMETER SCOPES						
KIT-0303-120	Purple	3.0 mm	2.0-2.2 mm	120 cm	Bronch/Cystos	50
KIT-0303-230	Purple	3.0 mm	2.0-2.2 mm	230 cm	Slim/Ped. Gastros	50
GI SCOPES						
KIT-0505-230	Blue	5.0 mm	2.0-4.5 mm	230 cm	GI Scopes	50
KIT-0505-250	Blue	5.0 mm	> 2.8 mm	250 cm	Enteroscopes	50


non-sterile

Note - all PMC Brush Kits include VCB-0511-15 5/11mm valve & port cleaning brush.

FLEXIBLE ENDOSCOPE BRUSHES

Single use

Single-ended Channel Brushes



CODE	SHEATH COLOR	BRUSH DIAMETER	COMPATIBLE CHANNEL	LENGTH	SCOPE TYPE	PKG QTY
CB-0002-100	White	2.0 mm	1.2-1.5 mm	100 cm	Microchannel	50
CB-0003-120	Purple	3.0 mm	2.0-2.2 mm	120 cm	Bronch/Cystos	50
CB-0005-230	Blue	5.0 mm	2.0-4.5 mm	230 cm	GI Scopes	50

Non-sterile

Double-ended Channel Brushes






CODE	DESCRIPTION	BRUSH DIAMETER	COMPATIBLE CHANNEL	LENGTH	PKG QTY
CB1015	Yellow sheath	2.0 mm	1.2-1.5 mm	100 cm	100
4665	Green sheath	5.0-6.0 mm	> 2.0 mm	100 cm	25
CB2024-L	White sheath	3.0 mm	2.0-2.4 mm	230 cm	50
CB2845	Blue sheath	5.0 mm	2.8-4.5 mm	230 cm	50
4660	Blue sheath	5.0-6.0 mm	> 2.0 mm	260 cm	25

Non-sterile

FLEXIBLE ENDOSCOPE BRUSHES

Single use

Pull Through Brush Models



Caterpillar Brush™
30 cm brush head length

CODE	DESCRIPTION	BRUSH DIAMETER	COMPATIBLE CHANNEL	LENGTH	PKG QTY
CAT-1	Blue sheath	-	2.7-4.5 mm	230 cm	100
CAT-2	Green sheath	-	2.0-2.2 mm	120 cm	100

Non-sterile



dragontail™ Brush
Microfiber strands

CODE	DESCRIPTION	BRUSH DIAMETER	COMPATIBLE CHANNEL	LENGTH	PKG QTY
DT2022	Teal sheath	-	2.0-2.2 mm	110 cm	100
DT2832	Purple sheath	-	2.8-3.2 mm	230 cm	100
DT3542	Pink sheath	-	3.5-4.2 mm	230 cm	100

Non-sterile

FLEXIBLE ENDOSCOPE BRUSHES

Reusable

Channel Brushes

CODE	BRUSH DIAMETER	COMPATIBLE CHANNEL	SHEATH LENGTH	PKG QTY
8202	0.8 mm	> 1.0 mm	120 cm	2
8204	2.0 mm	> 1.4 mm	120 cm	2
8206	3.0 mm	> 1.8 mm	120 cm	2

Non-sterile

FLEXIBLE ENDOSCOPE BRUSHES
Reusable

Valve Brushes



CODE	DESCRIPTION	BRUSH DIAMETER	BRUSH LENGTH	PKG QTY
8500	Red, double-ended	12.0 & 4.0 mm	2.5 & 1 cm	5
DV-82	Black, double-ended	11.0 & 4.0 mm	4 & 1 cm	5
SV-80	Black, single-ended	11.0 mm	4 cm	5
				Non-sterile

VALVE MANAGEMENT
Reusable

Mesh Bag

Durable, lightweight, polyester.
Compatible with high level disinfectants, gas plasma, and ETO sterilization.



CODE	DESCRIPTION	DIMENSIONS	PKG QTY
MB-0305	Lightweight mesh with drawstring closure	3″ x 5″	10
			Non-sterile

BEDSIDE PRE-CLEAN SUPPLIES

Single use

STEP 1:



Use the Krinkle Pad to wipe down the insertion tube.



STEP 2:



Suction the remaining fluid directly from the pack to clear out the channel.

Immediate bedside cleaning

First Step Bedside Kits simplify the endoscope cleaning process and prevent bioburden from drying and solidifying in scope channels ensuring CSA pre-cleaning standards are met.

Convenient – increases staff compliance & turnover time.

Ready to use – properly mixed formula every time.

Safe for use on all flexible scopes.

User-friendly, space efficient pouches are easily transportable.

Available in Ready to Use & Add Water Kits.

CODE	DESCRIPTION	VOLUME OF SOLUTION	PKG QTY
EP-3	READY TO USE Bedside Pre-Clean Kit with Krinkle Pad	100 ml	100
EP-4	READY TO USE Bedside Pre-Clean Kit with Krinkle Pad	200 ml	70
EP-4-30	READY TO USE Bedside Pre-Clean Kit with Krinkle Pad	200 ml	Pallet = 2100
EP-6	READY TO USE Bedside Pre-Clean Kit with Krinkle Pad	500 ml	30
EP-6W	ADD WATER Bedside Pre-Clean Kit with Krinkle pad	500 ml	120

Non-sterile

SCOPE CLEANING & TIP PROTECTION

Single use




Enzymatic Detergent Cleaning Kit

Pre-measured concentrate for proper dilution
(4 & 6 Gallons Kits).
1 package per scope – no need to reach for detergent pump.
Krinkle pad included.

CODE	DESCRIPTION	PKG QTY
EDK04	Simple2 Enzymatic Detergent, mixes with 4 gallons of water, one Krinkle Hand Pad	100
EP15G	Simple2 Enzymatic Detergent, mixes with 6 gallons of water, one Krinkle Hand Pad	100

Non-sterile



Endo Boot Cover

Star-shaped design opening for custom fit.
Protects distal end of scopes during storage or transport.



Non-Absorbent Scope Bumper

Does not capture or hold residual moisture that may lead to the formation of bacteria colonies.

GI Labs have reported up to a 50% decrease in bending rubber replacement with tip protection.

CODE	DESCRIPTION	PKG QTY
BC-50	Endo-Boot Cover - Large	50
BCB- 50	Endo-Boot Cover - Small	50
NA2500	Non-Absorbent Scope Bumper	200

Non-sterile


WIPES & CLEANING PADS

Single use




Clear Path™ Microfiber Wipes

Capture spores & other viable pathogens at a microscopic level.
Non-woven, lint-free, low-particulate split microfiber technology.
Compatible with any brand of surface disinfectant.
Proven to remove superbug pathogens C. Diff., VRE, CRE & Norovirus.



Draco™ Deep Cleaning Pads

Scrubbing Foam Side – extreme cleaning at a microscopic level.
Removes stubborn biofilm and adhered contaminants.
Split Microfiber Side – captures microscopic particles as small as 4 microns.
Inner Contour Foam Layer – wraps around instruments with ease.
Non-abrasive and safe for use on delicate optics.



Dual-Textured Enzymatic Cleaning Pads

Rough and smooth side accommodate all instruments.
Wet with tap water to activate detergent.
Simple2 Multi-tiered Enzymatic Detergent effectively dissolves blood, fat, tissue, protein, and most other forms of organic material.
Non-toxic, latex-free, neutral pH and safe for all surgical instruments.



Endo Scrub

Effective manual cleaning.
Easy to hold.



Dry Krinkle Pads

Remove debris with ease.
Non-abrasive.

CODE	DESCRIPTION	SIZE	PKG QTY
W5412-30	Microfiber Wipes (20 waterproof bags of 30 wipes)	8" x 12"	600
W5412-150	Microfiber Wipes (8 bags of 150 wipes)	10" x 12"	1200
HY0305	DRACO Deep Cleaning Pads	3.5" x 5"	100
SC-50	Endo Scrub – Dry Pad	9.5 cm L	50
3550DT	Dual Textured Enzymatic Cleaning Pads	3.5" x 5"	100
3550K	Dry Krinkle Cleaning Pad	3.5" x 5.5"	200

Non-sterile

SCOPE TRANSPORT CONTAINERS

Single use



OT1230

Features:

Color-coded reversible lids – CLEAN/SOILED.

Rigid construction protects scopes.

Biodegradable - made of a plant-based material.

Space Efficient – Stackable & Nestable.

No reprocessing of bins – No cross contamination.

3 Sizes to fit all scope sizes – small, standard & plus.



OT1016



OP1330

Made of a 100% Renewable Resource

CODE	DESCRIPTION	DIMENSIONS L x W	PKG QTY
OT1016	Oasis® Small Scope Transport Tray – Reversible lid No Reservoir	16" x 10"	50
OT1230NH	Oasis® Scope Transport Tray - Reversible lid Built in Reservoir	20.5" x 16.5"	30
OT1230NH-16	Oasis® Scope Transport Tray - Reversible lid Built in Reservoir - 16 Box Pallet	20.5" x 16.5"	480
OP1330	Oasis® PLUS Scope Transport Tray – Reversible lid No Reservoir	21.5" x 19.5"	30
OP1330-16	Oasis® PLUS Scope Transport Tray – Reversible lid No Reservoir - 16 Box Pallet	21.5" x 19.5"	480

Non-sterile

TRANSPORT CARTS

Single use



Dispensing Dolly



Collection Unit



TC06-R



TC04-G

OASIS® carts:

Color coded carts for clean **GREEN** and soiled **RED** scopes.
Green Dispensing Dolly holds 50 Clean trays & lids.
Red Collection unit holds soiled trays safely in biodegradable liners.
All feature non-marking swivel caster wheels & stainless steel base.

CODE	DESCRIPTION	DIMENSIONS W x D x H	COLOUR	PKG QTY
DD40-G	Oasis® & Oasis® PLUS Dispensing Dolly Removable tabletop, 50 tray capacity	23" x 21" x 36"	Green	1
CU10-R	Oasis® & Oasis® PLUS Collection Unit	23" x 21" x 36"	Red	1
TC04-G	Oasis® & Oasis® PLUS Transport Cart, 4 tray capacity	23" x 21" x 32"	Green	1
TC04-R	Oasis® & Oasis® PLUS Transport Cart, 4 tray capacity	23" x 21" x 32"	Red	1
TC06-G	Oasis® & Oasis® PLUS Transport Cart, 6 tray capacity	23" x 21" x 44"	Green	1
TC06-R	Oasis® & Oasis® PLUS Transport Cart, 6 tray capacity	23" x 21" x 44"	Red	1
TC08-G	Oasis® & Oasis® PLUS Transport Cart, 6 tray capacity	23" x 21" x 56"	Green	1
TC08-R	Oasis® & Oasis® PLUS Transport Cart, 8 tray capacity	23" x 21" x 56"	Red	1
TP-50	Oasis® & Oasis® PLUS Transport Top Cover	23"W x 21"D		1
4253-100	Biodegradable clear plastic liners for the Oasis® Tray Collection Unit	45 Gallon		100

STARTER KIT Includes 360 trays & rack – Contact for Quote & Order Code

TRANSPORT TRAYS

Reusable



Durable, strong, light-weight, plastic polymer.
Lid locks down to minimize spills or splashing.
Secured tray protects hands and eyes.
Removable insert tray for soaking and safe transferring of contaminated instruments.
Clearly identified with biohazard symbols and colours.

CODE	DESCRIPTION	DIMENSIONS L x W x H
8600	Tray with Instrument Soaking Insert	12.6 x 23.5 x 5.8"



Made of high temperature plastic.
High durability against heat and chemicals.
Transparent amber colour.
Easy locking, removable lid with removable silicone lining.

CODE	DESCRIPTION	DIMENSIONS L x W x H
IN-PIC-11	Large Tall	20.6 x 15.6 x 7.5"
IN-PIC-12	Large Short	20.6 x 15.6 x 5.3"
IN-PIC-13	Small Tall	13.6 x 15.6 x 7.5"
IN-PIC-14	Small Short	13.6 x 15.6 x 4.3"
IN-PIC-15	Transport container, XL short	20.5 x 12.625 x 3.75"
IN-PIC-16	Transport container, XL tall	20 x 12.25 x 6.125"

TRANSPORT BAGS

Single use



SIGNAL ™ Transport Bags

CLEARLY IDENTIFIES soiled scopes, protecting staff & patients.

COMPOSTABLE (ASTM D6400 standards).

TWO SIZES available.

PRE-APPLIED LABEL option:
Place HLD instrument inside bag & seal.
Record pkg. prep date & preparer's name on pre-applied HLD label.
After procedure, remove HLD label to expose Biohazard legend.
Place contaminated instrument in bag & reseal for transport.







Biodegradable & Standard Soiled Scope Bags

Kangaroo Pad

Cinch table cover into soiled scope containment bag.

CODE	DESCRIPTION	DIMENSIONS W x L	PKG QTY
KP100	Kangaroo Pad – Blue	41" x 36"	100
SB24W	Soiled Scope Bag – White	24" x 24"	100
HLD22	Biodegradable Soiled Scope Bag – Green	22" x 22"	100
ITB1218	SIGNAL™ – Compostable Transport Bag	12" x 18"	200
ITB2022	SIGNAL™ – Compostable Transport Bag	20" x 22"	200
ITB1218-H	SIGNAL™ – Compostable Transport Bag w/ HLD label	12" x 18"	200
ITB2022-H	SIGNAL™ – Compostable Transport Bag w/ HLD label	20" x 22"	200

Non-sterile

COUNTDOWN™ SCOPE TAGS

Single use





MONTH/DAY or DAILY Tags

Information you need, at a glance.

Clearly labeled scope tags quickly identify when scopes were reprocessed, aiding in compliance to industry standards for storage intervals.

Convenient pressure sensitive adhesive strip.
Color-coded day of week.
Latex-free.

BEST PRACTICE

CSA Standard Z314-18 Canadian Medical Device Reprocessing indicates in section 12.5.4 Endoscope shelf life 1.3.1.3 There shall be a process in place to identify the date of reprocessing for each endoscope.

CODE	DESCRIPTION	DIMENSIONS	PKG QTY
CD300	Month/Day – 2 pks of 100	12" x 1.75"	200
CDM610	MONDAY – 4 pks of 100	12" x 1.75"	400
CDT620	TUESDAY – 4 pks of 100	12" x 1.75"	400
CDW630	WEDNESDAY – 4 pks of 100	12" x 1.75"	400
CDTH640	THURSDAY – 4 pks of 100	12" x 1.75"	400
CDF650	FRIDAY – 4 pks of 100	12" x 1.75"	400
CDS660	SATURDAY – 4 pks of 100	12" x 1.75"	400
CDSUN670	SUNDAY– 4 pks of 100	12" x 1.75"	400
CDCOMB700	Each Weekday – 100 tags	12" x 1.75"	700

GI DEVICES & ACCESSORIES

POLYP TRAPS	24
BITE BLOCKS	25
CAPSULE ENDOSCOPY	26
RETRIEVAL DEVICES	27
RETRIEVAL FORCEPS - SPECIALTY	28
POLYPECTOMY SNARES	29
INJECTION NEEDLES	30
ENDOSCOPIC HEMOCLIP	31
CYTOLOGY & MICROBIOLOGY BRUSHES	33
BALLOON DILATION SYSTEM - GI	34
BALLOON DILATORS	35
ACHALASIA DILATORS	37
BILIARY DILATORS	38
BILIARY STENTS	39
PANCREATIC FLEXI-STENTS	40
CATHETERS	42

Keir Surgical Endoscope Care

KS COMPLETE CARE

Do complicated Endoscope Repairs cause you headaches?
Let us simplify that process.
We are pleased to offer flexible endoscope repairs, all in one convenient program.

OEM CORRECT™



WHAT WE OFFER:

- **Fast turnaround:** Repairs completed at our state of the art facility in Vaughan, ON
- **High Quality:** ISO 13485 certified & our team of technicians have over 250 years of combined OEM Correct™ endoscope repair experience!
- **Cost savings:** Significant savings compared to OEM repairs.

Comprehensive **Save Time & Money** **Trusted**

For more information, contact us at:
705.201.1188 | repairs@keirsurgical.com

POLYP TRAPS

Single use

Accommodate a single polyp or many, save budget dollars by not discarding empty chambers. Stock a variety!

Single Chamber

Screw cap easily removes when suction engaged.
With or without blue tip, accommodating various endoscope suction connections.
Multi-use trap can also collect ENT sinus polyps, arthroscopy shavings & urology stones.



CODE	DESCRIPTION	PKG QTY
PT-25	Blue Connector Tip	25
PT-25X	NO Blue Tip	25
PT-100	Blue Connector Tip	100
PT-100X	NO Blue Tip	100

Non-sterile

Multiple Chambers



CODE	DESCRIPTION	PKG QTY
VDK-PT-05	5 Chambers	10

Non-sterile

BITE BLOCKS

Single use









CODE	DESCRIPTION	INNER DIAMETER	PKG QTY
VDK-KD-B1	Green with white pliable strap	20 mm/60 FR	100
EN-210	Pink with blue pliable strap	16 mm/48 FR	25
EN-300	Orange strapless	16 mm/48 FR	25
EN-310	Blue with blue pliable strap	21 mm/63 FR	25
EN-320	Purple with velcro strap	21 mm/63 FR	25

Non-sterile

CAPSULE ENDOSCOPY

Single Use

CapsoCam Plus®



Desiged to improve diagnostic capability by providing a full 360º panoramic direct lateral view of the small bowel mucosa.

No external data recorder

- Increased patient comfort.
- Safe to use with patients with pacemakers or other electronic implants.
- No risk of data loss with patients with high BMI.
- Ideal for remote patients.

Full 360° View



CapsoCam Plus image of the small bowel mucosa viewed at 1152 x 212 pixels.





Competitor *CapsoCam Plus*

CODE	DESCRIPTION	PKG QTY
SV	CAPSOCAM PLUS CAPSULE ENDOSCOPE & RETRIEVAL KIT, EACH	1
CDAS	CAPSOCAM PLUS CAPSULE DATA ACCESS	1
		Sterile

RETRIEVAL DEVICES

Reusable

Pronged Retrievers



**Replacement handles available*

CODE	DESCRIPTION	OPEN DIAM.	SHEATH DIAM.	LENGTH OVERALL	PKG QTY
5004	3 Prongs	20 mm	1.8 mm	150 cm	1

Sterile

Basket Retrievers



**Replacement handles available*

CODE	DESCRIPTION	OPEN DIAM.	SHEATH DIAM.	LENGTH OVERALL	PKG QTY
4816	3 Wires	20 mm	2.3 mm	220 cm	1
4814	4 Wires	15 mm	1.8 mm	150 cm	1
4832	4 Wires	20 mm	2.3 mm	220 cm	1
4822	4 Wires	15 mm	1.8 mm	150 cm	1

Sterile

RETRIEVAL FORCEPS - SPECIALTY

Single Use

Rotatable Net Retriever



CODE	NET SIZE	SHEATH DIAM.	LENGTH	PKG QTY
VDK-FG-25-230-C3	30 x 60 mm	2.5	230 cm	10
VDK-FG-25-230-C4	40 x 75 mm	2.5	230 cm	10

Sterile

Looped Retriever

Unique spring-loaded handle



CODE	DESCRIPTION	OPEN DIAM.	MIN. CHANNEL	LENGTH	PKG QTY
5030S	3 Looped Prongs	2.5 cm	≥ 2.8 mm	220 cm	5

Sterile

Grasping Forceps

Diversatek™ EndoGrip




CODE	DESCRIPTION	SHEATH DIAM.	MIN. CHANNEL	LENGTH	PKG QTY
1177-01	Rat Tooth Grasper with Alligator Jaws	2.3 mm	≥ 2.8 mm	230 cm	5

Sterile

*Additional styles may be available. Please contact your representative for details.

POLYPECTOMY SNARES

Single-Use - Use with Olympus Active Cords



CODE	DESCRIPTION	LOOP SIZE	SHEATH DIAM.	LENGTH	PKG QTY
4561-OLY	Hexagonal	4.0 x 3.0 cm	2.3 mm	220 cm	5
4562-OLY	Crescent	5.0 x 2.5 cm	2.3 mm	220 cm	5
4563-OLY	Standard Oval	5.0 x 2.5 cm	2.3 mm	220 cm	5
4564-OLY	Maxi-Oval	5.5 x 3.5 cm	2.3 mm	220 cm	5
4565-OLY	Mini-Oval	3.5 x 2.0 cm	2.3 mm	220 cm	5
4567-OLY	Micro-Oval	2.5 x 1.5 cm	2.3 mm	220 cm	5

Sterile

Short Throw Handle




CODE	DESCRIPTION	LOOP SIZE	SHEATH DIAM.	LENGTH	PKG QTY
4573-OLY	Standard Oval	5.0 x 2.5 cm	2.3 mm	220 cm	5

Sterile

INJECTION NEEDLES

Single-Use



Basic operation for ease of use in emergency situations.

Tip guard and needle lock prevent potential punctures to scope channel.


CODE	DESCRIPTION	SHEATH DIAMETER	LENGTH	PKG QTY
GASTROSCOPY LENGTH				
4702	23 gauge/ 5.0 mm	1.8 mm	200 cm	1
4704	25 gauge/ 5.0 mm	1.8 mm	200 cm	1
4706	23 gauge/ 5.0 mm	2.3 mm	200 cm	1
4708	25 gauge/ 5.0 mm	2.3 mm	200 cm	1
4714	25 gauge/ 4.0 mm	1.8 mm	200 cm	1
4718	25 gauge/ 4.0 mm	2.3 mm	200 cm	1
COLON/ENTEROSCOPY LENGTH				
4720	23 gauge/ 5.0 mm	2.3 mm	260 cm	1
4721	23 gauge/ 4.0 mm	2.3 mm	260 cm	1
4722	25 gauge/ 5.0 mm	2.3 mm	260 cm	1
4723	25 gauge/ 4.0 mm	2.3 mm	260 cm	1

Sterile

ENDOSCOPIC HEMOCLIP

Single-Use

One-to-one rotation
Accurate positioning
Handle safety lock design
Repositionable clips






CODE	CLIP OPENING SIZE	WORKING LENGTH	ENDOSCOPE CHANNEL	PKG QTY
VDK-CD-11-230-O (Uncoated)	11 mm	230 cm	≥ 2.8 mm	5
VDK-CD-11-230-P	11 mm	230 cm	≥ 2.8 mm	5
VDK-CD-13-230-P	13 mm	230 cm	≥ 2.8 mm	5
VDK-CD-16-230-P	16 mm	230 cm	≥ 2.8 mm	5
				Sterile

BIOPSY FORCEPS

Single-use

4-bar linkage design, offering excellent clamping strength.
Smooth jaws, decreasing the friction between the working channel.
Standard cup design with proper sampling capacity and minimal tissue harm.
Many different patterns available upon request!



STRONG BITE FORCE
Easily grasp the mucosal tissue, reducing risks of tearing and bleeding.



CODE	JAW STYLE	JAW DIAM.	WORKING CHANNEL	LENGTH	PKG QTY
VDK-FB-18-105-O-P-A1	Without spike, coated, with alligator teeth	1.8 mm	≥ 2.0 mm	105 cm	10
VDK-FB-18-105-O-P-B1	Without spike, coated	1.8 mm	≥ 2.0 mm	105 cm	10
VDK-FB-18-105-W-P-B1	With spike, coated	1.8 mm	≥ 2.0 mm	105 cm	10
VDK-FB-18-160-O-P-B1	Without spike, coated	1.8 mm	≥ 2.0 mm	160 cm	10
VDK-FB-23-160-W-O-A1	With spike, uncoated, with alligator teeth	2.3 mm	≥ 2.8 mm	160 cm	10
VDK-FB-23-230-O-P-A1	Without spike, coated, with alligator teeth	2.3 mm	≥ 2.8 mm	230 cm	10
VDK-FB-23-230-O-P-B1	Without spike, coated	2.3 mm	≥ 2.8 mm	230 cm	10
VDK-FB-23-230-W-P-A1	With spike, coated, with alligator teeth	2.3 mm	≥ 2.8 mm	230 cm	10
VDK-FB-23-230-W-P-B1	With spike, coated	2.3 mm	≥ 2.8 mm	230 cm	10


Sterile

CYTOLOGY & MICROBIOLOGY BRUSHES

Single-Use

Cytology Brushes

Simple handle enables easy rotation.
Brush designed for maximum yield with minimum deposits left in sheath.
Radiopaque ball tip prevents mucosal trauma.
High tempered stainless-steel wire prevents kinks.




CODE	BRUSH DIAMETER	SHEATH DIAMETER	LENGTH	PKG QTY
4202	1.0 mm	1.7 mm	120 cm	10
4204	2.0 mm	1.7 mm	120 cm	10
4206	3.0 mm	1.7 mm	120 cm	10
4102	3.0 mm	1.7 mm	160 cm	10
4104	3.0 mm	1.7 mm	200 cm	10
4106	3.0 mm	1.7 mm	230 cm	10
4108	3.0 mm	2.3 mm	160 cm	10
4110	3.0 mm	2.3 mm	200 cm	10
4112	3.0 mm	2.3 mm	260 cm	10

Sterile

Microbiology Brushes

Design prevents brush contamination.
Spring loaded handle for automatic seal of distal tip.
Handle locking device for open brush position.
Brush configuration increases specimen yield.




CODE	BRUSH DIAMETER	SHEATH DIAMETER	LENGTH	PKG QTY
4352	2.0 mm	1.7 mm	220 cm	5
4310	1.0 mm	1.7 mm	100 cm	5

Sterile

BALLOON DILATION SYSTEM - GI

Uniform dilation achieved with gradual inflation pressure.
All parts sold individually and compatible with other dilation systems.

SIMPLE 2-PART SETUP!




1. Balloon Dilators



Stylet and Over-the-Wire balloon dilators.
See following pages for selection.

2. Inflation Monitor



Luminescent gauge.
Universal fittings adapt to other dilators.

CODE	PKG QTY
4630	1

Non-Sterile

BALLOON DILATORS

Stylet Wire - Single Use



Syringe and monitor not included

Uniform dilation achieved with gradual inflation pressure.

Compatible with most flexible GI scopes.


Interchangeable with other dilator systems.

CODE	INFLATED OUTER DIAMETER	BALLOON LENGTH	SHEATH DIAMETER	WORKING LENGTH	PKG QTY
3710	14 mm / 42 FR	8 cm	2.3 mm	240 cm	1
3712	16 mm / 48 FR	8 cm	2.3 mm	240 cm	1
3714	18 mm / 54 FR	8 cm	2.3 mm	240 cm	1
3716	20 mm / 60 FR	8 cm	2.3 mm	240 cm	1
3732	8 mm / 24 FR	5.5 cm	2.3 mm	240 cm	1
3734	10 mm / 30 FR	5.5 cm	2.3 mm	240 cm	1
3736	12 mm / 36 FR	5.5 cm	2.3 mm	240 cm	1
3738	14 mm / 42 FR	5.5 cm	2.3 mm	240 cm	1
3740	16 mm / 48 FR	5.5 cm	2.3 mm	240 cm	1
3742	18 mm / 54 FR	5.5 cm	2.3 mm	240 cm	1
3744	20 mm / 60 FR	5.5 cm	2.3 mm	240 cm	1
3752	8 mm / 24 FR	3 cm	2.3 mm	240 cm	1
3754	10 mm / 30 FR	3 cm	2.3 mm	240 cm	1
3756	12 mm / 36 FR	3 cm	2.3 mm	240 cm	1
3758	14 mm / 42 FR	3 cm	2.3 mm	240 cm	1
3760	16 mm / 48 FR	3 cm	2.3 mm	240 cm	1
3762	18 mm / 54 FR	3 cm	2.3 mm	240 cm	1
3764	20 mm / 60 FR	3 cm	2.3 mm	240 cm	1

Non-Sterile

BALLOON DILATORS

Over-the-Wire - Single Use



Atraumatic, Teflon-coated guide wire passes easily through working channels and gradually penetrates narrow strictures.

Sizes for standard and therapeutic scopes.

Interchangeable with other dilation systems.

Visible under fluoroscopy.

CODE	INFLATED OUTER DIAMETER	BIOPSY CHANNEL	BALLOON LENGTH	SHEATH DIAMETER	WORKING LENGTH	PKG QTY
3152	8 mm / 24 FR	2.8 mm	3 cm	2.3 mm	240 cm	1
3156	10 mm / 30 FR	2.8 mm	3 cm	2.3 mm	240 cm	1
3158	12 mm / 36 FR	2.8 mm	3 cm	2.3 mm	240 cm	1
3206	14 mm / 42 FR	3.7 mm	3 cm	2.3 mm	240 cm	1
3208	16 mm / 48 FR	3.7 mm	3 cm	2.3 mm	240 cm	1
3210	18 mm / 54 FR	3.7 mm	3 cm	2.3 mm	240 cm	1
3212	20 mm / 60FR	3.7 mm	3 cm	2.3 mm	240 cm	1
3102	6 mm / 18 FR	2.8 mm	4 cm	2.3 mm	240 cm	1
3104	8 mm / 24 FR	2.8 mm	4 cm	2.3 mm	240 cm	1
3106	10 mm / 30 FR	2.8 mm	4 cm	2.3 mm	240 cm	1
3202	16 mm / 48 FR	3.7 mm	4 cm	2.3 mm	240 cm	1
3204	20 mm / 60 FR	3.7 mm	4 cm	2.3 mm	240 cm	1
3250	8 mm / 24 FR	3.7 mm	8 cm	2.3 mm	240 cm	1
3254	10 mm / 30 FR	3.7 mm	8 cm	2.3 mm	240 cm	1
3258	12 mm / 36 FR	3.7 mm	8 cm	2.3 mm	240 cm	1
3260	14 mm / 42 FR	3.7 mm	8 cm	2.3 mm	240 cm	1
3262	16 mm / 48 FR	3.7 mm	8 cm	2.3 mm	240 cm	1
3264	18 mm / 54 FR	3.7 mm	8 cm	2.3 mm	240 cm	1
3266	20 mm / 60 FR	3.7 mm	8 cm	2.3 mm	240 cm	1

Sterile

Guide Wires

Sold separately

CODE	DESCRIPTION	LENGTH	OUTER DIAMETER	PKG QTY
3412	Teflon-coated, 7 cm floppy tip	400 cm	.035”	1

Sterile

ACHALASIA DILATORS

Single Use



Syringe and monitor not included
Use with .038 Guide Wire

Uniform shaft diameter provides greater rigidity when passed through the lower esophageal sphincter.

Sturdy, specially treated polyethylene promotes durability.

Visible under fluoroscopy.

CODE	INFLATED OUTER DIAMETER	BALLOON LENGTH	SHEATH DIAMETER	WORKING LENGTH	PKG QTY
3302	30 mm / 90 FR	10 cm	5.3 mm	100 cm	1
3304	35 mm / 105 FR	10 cm	5.3 mm	100 cm	1
3306	40 mm / 120 FR	10 cm	2.3 mm	100 cm	1

Sterile

Guide Wires

Sold separately

CODE	DESCRIPTION	LENGTH	OUTER DIAMETER	PKG QTY
3406	Teflon-coated, 7 cm floppy tip	260 cm	.038"	1

Sterile

BILIARY DILATORS

Single Use



Sturdy, specially treated polyethylene promotes successful dilation.

Teflon-coated guide wire passes easily through working channel.

Interchangeable with other dilation systems.

Visible under fluoroscopy.

Over-the-Wire - Dilator

CODE	INFLATED OUTER DIAMETER	BALLOON LENGTH	SHEATH DIAMETER	WORKING LENGTH	BIOPSY CHANNEL	PKG QTY
3180	4 mm / 12 FR	2 cm	2.3 mm	180 cm	2.8 mm	1
3182	6 mm / 18 FR	2 cm	2.3 mm	180 cm	1.8 mm	1
3184	8 mm / 24 FR	3 cm	2.3 mm	180 cm	2.8 mm	1

Sterile

Guide Wires


Sold separately

CODE	DESCRIPTION	LENGTH	OUTER DIAMETER	PKG QTY
3412	Teflon-coated, 7 cm floppy tip	400 cm	.035"	1
3413	Teflon-coated, 7 cm floppy tip	480 cm	.035"	1

Sterile

BILIARY STENTS

Single Use



Also available in kits with guide wire and push catheter

Unique coating for enhanced biocompatibility.

Extra flexible material reduces risk of duct trauma.

MR safe.

Pigtail stent – Equal radius on ends.


Six equally spaced side ports for superior drainage.

Amsterdam Stent – internal/external retention flanges for easy positioning.

CODE	DESCRIPTION	DIAMETER	LENGTH	WORKING LENGTH	PKG QTY
6102	Pigtail Stent	7 FR	4 cm	240 cm	1
6103	Pigtail Stent	7 FR	7 cm	240 cm	1
6104	Pigtail Stent	7 FR	10 cm	240 cm	1
6106	Pigtail Stent	10 FR	4 cm	240 cm	1
6107	Pigtail Stent	10 FR	7 cm	240 cm	1
					Sterile

PANCREATIC FLEXI-STENTS

Single use



Also available in kits with push catheter

Freeman-Aliperti Straight Style

Designed to bridge sphincter into distal duct.

Double external flanges and distal anchoring flange prevent stent migration.

Pliable material & softer flanges minimize duct trauma & promote patient comfort.

External marks for proper placement.

Compatible with current exchange systems.

MR Safe.


CODE	STENT DIAMETER	STENT LENGTH	GUIDE WIRE COMPATIBILITY	PKG QTY
6341	4 FR	2 cm	.025"	1
6342	4 FR	3 cm	.025"	1
6351	5 FR	2 cm	.025"	1
6352	5 FR	3 cm	.025"	1
				Sterile

Push/Guide Catheters – Pancreatic & Biliary Stents

CODE	DESCRIPTION	MINIMUM CHANNEL DIAM.	LENGTH	PKG QTY
6138	Push Catheter for 3 FR stents	2.8 mm	170 cm	1
6136	Push Catheter for 4 & 5 FR stents	2.8 mm	170 cm	1
6132	Push Catheter for 7 FR stents	2.8 mm	170 cm	1
6134	Push Catheter for 10 FR stents	2.8 mm	170 cm	1
				Sterile

PANCREATIC FLEXI-STENTS

Single use



Also available in kits with push catheter

Freeman-Aliperti Straight Style

Pliable material easily follows contours of tortuous ducts – less trauma, lower risk of post ERCP pancreatitis.
Softer retention flanges minimize duct trauma.
Small Pigtail radius prevents interference with duodenal peristalsis.
Compatible with current exchange systems.
MR Safe.

CODE	DESCRIPTION	DIAMETER	LENGTH	GUIDE WIRE COMPATIBILITY	PKG QTY
6532	Without Flange	3 FR	5 cm	.021″	1
6541	With Flange	4 FR	3 cm	.025″	1
6542	With Flange	4 FR	5 cm	.025″	1
6544	Without Flange	4 FR	7 cm	.025″	1
6545	Without Flange	4 FR	9 cm	.025″	1
6546	Without Flange	4 FR	11 cm	.025″	1
6551	With Flange	5 FR	3 cm	.035″	1
6552	With Flange	5 FR	5 cm	.035″	1
6554	With Flange	5 FR	7 cm	.035″	1
6555	With Flange	5 FR	9 cm	.035″	1
6561	With Flange	5 FR	3 cm	.035″	1
6562	Without Flange	5 FR	5 cm	.035″	1
6564	Without Flange	5 FR	7 cm	.035″	1
6565	Without Flange	5 FR	9 cm	.035″	1
6566	Without Flange	5 FR	12 cm	.035″	1
6567	Without Flange	5 FR	15 cm	.035″	1
6571	With Flange	7 FR	3 cm	.035″	1
6572	With Flange	7 FR	5 cm	.035″	1
6574	With Flange	7 FR	7 cm	.035″	1
6575	With Flange	7 FR	9 cm	.035″	1
6576	With Flange	7 FR	11 cm	.035″	1
6582	Without Flange	7 FR	5 cm	.035″	1
6584	Without Flange	7 FR	7 cm	.035″	1
6585	Without Flange	7 FR	9 cm	.035″	1

Sterile

CATHETERS

Single use



Spray Catheters

Precision spray nozzle delivers a powerful 360-degree spray mist pattern for use in Chromoendoscopy, as well as for washing and power irrigation.

Consistent 360-degree spray mist pattern.

Compatible with various irrigators and Hobbs Vacu-Lock syringe.

CODE	SHEATH DIAMETER	LENGTH	PKG QTY
2190	2.6 mm	260 cm	5
			Sterile



Aspiration Catheters

Contaminant free collection of fluids.

Permanent or removable Luer-lock adapter.

Removable adapter allows scope removal during aspiration procedure.

Distal tip spiral hole pattern improves fluid collection.

PTFE catheter for ease of passage through channel.

Can be used with Hobbs Vacu-Lok Aspirating Syringe.

CODE	SHEATH DIAMETER	LENGTH	LUER-LOCK ADAPTER	PKG QTY
2181	2.3 mm	240 cm	Permanent	5
2182	2.3 mm	240 cm	Removable	5
				Sterile

PRODUCT INDEX

1177-01	28	**3744**	35	**5030S**	28	**ASP2**	7
2181	42	**3752**	35	**6102**	39	**BC-50**	15
2182	42	**3754**	35	**6103**	39	**BCB- 50**	15
2190	42	**3756**	35	**6104**	39	**CAT-1**	12
3102	36	**3758**	35	**6106**	39	**CAT-2**	12
3104	36	**3760**	35	**6107**	39	**CB1-110**	6
3106	36	**3762**	35	**6132**	40	**CB2-60**	6
3152	36	**3764**	35	**6134**	40	**CB-0002-100**	11
3156	36	**4102**	33	**6136**	40	**CB2-110**	6
3158	36	**4104**	33	**6138**	40	**CB2-200**	6
3180	38	**4106**	33	**6341**	40	**CB-0003-120**	11
3182	38	**4108**	33	**6342**	40	**CB-0005-230**	11
3184	38	**4110**	33	**6351**	40	**CB1015**	11
3202	36	**4112**	33	**6352**	40	**CB2024-L**	11
3204	36	**4202**	33	**6532**	41	**CB2845**	11
3206	36	**4204**	33	**6541**	41	**CD300**	21
3208	36	**4206**	33	**6542**	41	**CDAS**	26
3210	36	**4253-100**	18	**6544**	41	**CDCOMB700**	21
3212	36	**4310**	33	**6545**	41	**CDF650**	21
3250	36	**4352**	33	**6546**	41	**CD-HDMI**	6
3254	36	**4561-OLY**	29	**6551**	41	**CDM610**	21
3258	36	**4562-OLY**	29	**6552**	41	**CDS660**	21
3260	36	**4563-OLY**	29	**6554**	41	**CDSUN670**	21
3262	36	**4564-OLY**	29	**6555**	41	**CDT620**	21
3264	36	**4565-OLY**	29	**6561**	41	**CDTH640**	21
3266	36	**4567-OLY**	29	**6562**	41	**CD-USB**	6
3302	37	**4573-OLY**	29	**6564**	41	**CDW630**	21
3304	37	**4630**	34	**6565**	41	**CS9179**	8
3306	37	**4660**	11	**6566**	41	**CU10-R**	18
3406	37	**4665**	11	**6567**	41	**DD40-G**	18
3412	36	**4702**	30	**6571**	41	**DF100X**	7
3412	38	**4704**	30	**6572**	41	**DV-82**	13
3413	38	**4706**	30	**6574**	41	**EDK04**	15
3550DT	16	**4708**	30	**6575**	41	**EN-210**	25
3550K	16	**4714**	30	**6576**	41	**EN-300**	25
3710	35	**4718**	30	**6582**	41	**EN-310**	25
3712	35	**4720**	30	**6584**	41	**EN-320**	25
3714	35	**4721**	30	**6585**	41	**EP-3**	14
3716	35	**4722**	30	**8202**	12	**EP-4**	14
3732	35	**4723**	30	**8204**	12	**EP-4-30**	14
3734	35	**4814**	27	**8206**	12	**EP-6**	14
3736	35	**4816**	27	**8500**	13	**EP-6W**	14
3738	35	**4822**	27	**8600**	19	**EP15G**	15
3740	35	**4832**	27	**ASF2**	7	**FLT-Y30**	7
3742	35	**5004**	27	**ASO3**	7	**FLT-Y40ST**	7

PRODUCT INDEX

HLD22 20
HY0305 16
IN-PIC-11 19
IN-PIC-12 19
IN-PIC-13 19
IN-PIC-14 19
IN-PIC-15 19
IN-PIC-16 19
ITB1218 20
ITB1218-H 20
ITB2022 20
ITB2022-H 20
IV500 8
KIT-0002-100 10
KIT-0003-120 10
KIT-0005-230 10
KIT-0202-100 10
KIT-0303-120 10
KIT-0303-230 10
KIT-0505-230 10
KIT-0505-250 10
KP100 20
LL01-01 7
LL01-25 7
MB-0305 13
NA2500 15
OP1330 17
OP1330-16 17
OT1016 17
OT1230NH 17
OT1230NH-16 17
PT-25 24
PT-25X 24
PT-100 24
PT-100X 24
SB24W 20
SC-50 16
SDT01 7
SDT02 7
SV 26
SV-80 13
TC04-G 18
TC04-R 18
TC06-G 18
TC06-R 18
TC08-G 18
TC08-R 18
TP-50 18
VCB-10 10
VCB-0511-15 10
VDK-CD-11-230-O (UN-COATED) 31
VDK-CD-11-230-P 31
VDK-CD-13-230-P 31
VDK-CD-16-230-P 31
VDK-FB-18-105-O-P-A1 32
VDK-FB-18-105-O-P-B1 32
VDK-FB-18-105-W-P-B1 32
VDK-FB-18-160-O-P-B1 32
VDK-FB-23-160-W-O-A1 32
VDK-FB-23-230-O-P-A1 32
VDK-FB-23-230-O-P-B1 32
VDK-FB-23-230-W-P-A1 32
VDK-FB-23-230-W-P-B1 32
VDK-FG-25-230-C3 28
VDK-FG-25-230-C4 28
VDK-KD-B1 25
VDK-PT-05 24
VFR100 7
VFR100 7
W5412-30 16
W5412-150 16
WRC100 8
WRC500 8

NOTES

NOTES

KEIR SURGICAL

INSTRUMENTS FOR LIFE
DES INSTRUMENTS POUR LA VIE

Head Office
Tel: 800.663.4525 or 604.261.9596 | Fax: 604.261.9549
126-408 East Kent Avenue South, Vancouver, BC V5X 2X7

Eastern Regional Office
Tel: 800.663.4525 ou 514.277.5347 | Fascimile: 514.270.5340
9460 Route Transcanadienne, Saint-Laurent, QC H4S 1R7

info@keirsurgical.com | www.keirsurgical.com

PACIFIC SURGICAL

Keir Surgical is a Pacific Surgical Company
Keir Surgical est une filiale de Pacific Surgical